



Nymi Connect for Windows Release Notes

NCW 1.2.0

July 3, 2026
Version 1

Contents

Contents.....	2
Nymi Connect for Windows 1.2.0.....	3
New In This Release.....	4
Compatibility.....	4
Known Issues and Limitations	5

Table 1: Revision History

Version	Date	Revision History
1	July 3, 2026	First release of document

Nymi Connect for Windows 1.2.0

This release expands Nymi Connect for Windows (NCW) to support Java-based legacy applications, enabling passwordless “tap-and-go” authentication with the Nymi Band across a broader set of Manufacturing Execution Systems (MES).

This release enhances Nymi Connect for Windows (NCW) by adding support for Java-based legacy applications, significantly expanding the reach of Nymi’s passwordless authentication in pharmaceutical manufacturing environments

The NCW 1.2.0 is a minor release with Risk Level 0. It can be used for use in production deployments or GxP environments.

Nymi Connect for Windows components

NCW 1.2.0 is a single component release and includes the following component versions:

- Nymi Connect for Windows 1.2.0

Nymi Documentation Package

Documentation has been updated for NCW 1.2.0. For latest documentations visit the Nymi Support site.

- To access all NCW core documents, visit: [Nymi Connect for Windows Product Release Documentation](#)

NFC Reader Support

NFC readers were not explicitly tested under the present release. For a list of NFC readers supported in this release, please refer to the [Supported NFC Readers](#) document.

New In This Release

Nymi Connect for Java Based Applications

- Added support for Java-based MES applications.
- Extended NCW compatibility beyond native Windows and web apps to various Java UI frameworks.

Compatibility

Compatibility

NCW 1.2.0 is compatible with CWP 1.20.0 and above.

Known Issues and Limitations

Table 2: Known Issues in NCW 1.2.0

Jira ID	Summary	Justification
CONNECT-131	Auto-select JAB CanExecute probe stalls for seconds (~11s observed) on non-Java foreground windows (e.g. Edge), delaying/starving the UIA injector	Limited to the auto-select injection path when a non-Java application is in the foreground. May introduce intermittent injector selection delays, but functionality remains unaffected once the correct injector is selected. Low release risk.
CONNECT-136	uninstalling Java could not close Nymi Connect	Occurs only when Java is uninstalled while Nymi Connect is running. Impact is limited to an uncommon maintenance scenario and can be resolved by closing Nymi Connect before uninstalling Java.
NEM-3811	NCW does not display an error message when Auto-Fill fails due to JAB being disabled or unavailable.	Affects only cases where JAB is disabled after NCW has started. Restarting NCW refreshes the JAB status. Limited impact and low likelihood of occurrence.
NEM-3825	NCW: marshal band-message read as UTF-8 to avoid non-ASCII username corruption	Impacts only users with non-ASCII usernames in specific system encoding configurations. Password handling is unaffected and overall occurrence is expected to be low.
NEM-3832	Different Error Messages Are Displayed for Windows and Web Applications When Focus Changes During Auto-Fill	Limited to inconsistent error messaging when application focus changes during Autofill. Auto-Fill behaviour remains correct with no impact to functionality or security.

Table 3: Fixed Issues in NCW 1.2.0

Nil

Table 4: Known Issues in NCW 1.1.0

Issue Number	Description
NEM-3763	<p>[NCW] Tool tip showing incorrect error when NCW cannot reach Nymi Agent</p> <p>Workaround: When multiple issues occurs then oldest error is shown instead of newest error regarding network connectivity.</p>
NEM-3772	<p>In Edge, sometimes inject doesn't work, sometime it doesn't</p> <p>Workaround: 1) disable Autofill or 2) Do not move mouse from user name field and Autofill is popped-up</p>
NEM-3738	<p>[NCW] NES URL should not be case sensitive</p> <p>Workaround: for NES URL registry, using proper case sensitive nes service name will not cause any issue.</p>
NEM-3741	<p>[NCW] On init when NBE is not detected, notification says NCW is ready for use</p> <p>Workaround: when user hovers the mouse to NCW icon , User will see correct state of NCW.</p>
NEM-3750	<p>[NCW] User cannot zip up NCW client and NCS log folders</p> <p>Workaround: to collect logs manually selecting all files and place them to reachable location if needed.</p>



Nymi Inc. (Nymi) believes the information in this document is accurate as of its publication date. The information is subject to change without notice.

The information in this document is provided as-is and Nymi makes no representations or warranties of any kind. This document does not provide you with any legal rights to any intellectual property in any Nymi product. You may copy and use this document for your referential purposes.

This software or hardware is developed for general use in a variety of industries and Nymi assumes no liability as a result of their use or application. Nymi, Nymi Band, and other trademarks are the property of Nymi Inc. Other trademarks may be the property of their respective owners.

Published in Canada.

Nymi Inc.

Toronto, Ontario

www.nymi.com